A. Indications for Use

The Razor and Everest are indicated for use by dental technicians in the construction of custom made all ceramic restorations for anterior and posterior location.

Caution: Federal law restricts this device to sale by or on the order of dental health professionals.

B. Device Description

Razor and Everest are zirconia-based ceramic provided in various shapes such as round and square, and they are used to manufacture cores of all ceramic crowns, and are classified into ISO 6872 Type 2 Class 5. They are used to manufacture ceramic restorations through cutting process by dental MAD/MAM, computer-assisted design system, or CAD/CAM system. We offer 37 different shades A0, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, BL1, BL3, ML A1, ML A2, ML A3, ML A3.5, ML A4, ML B1, ML B2, ML B3, ML B4, ML C1, ML C2, ML C3, ML C4, ML D2, ML D3, ML D4, ML BL1, ML BL3 to meet the needs of different patients' tooth colors. It also offers various shape types to be used with various jigs for CAD/CAM milling machine.

C. Preparations

- (1) Look into the product to check if any damage, impurities, or malfunction exists.
- (2) Check if the exterior or shape of the product is appropriate for use.
- (3) Check that the product is the same as the one intended for use.

D. How to Use

- (1) Use the proper seating jig to fix the product to the equipment in use.
- If resin or wax is being used, bury it until the bridge or cap is completely filled up.
- If bond is being used to install in the jig, check the adhesiveness prior to use.
- (2) Operate the equipment and cut.
- The minimum thickness: before sintering, approx. 0.6 mm, and after sintering, approx. 0.5 mm.
- (3) After cutting, remove the bridge or cap carefully.
- (4) Use a drying machine or the like to dry the bridge or cap.
- (5) Put into the sintering furnace and bake.

- Recommended sintering schedule

Temperature	Time		
	Temperature/min	Running Time	Holding Time
*R.T. ~ 1000°C	17°C/min	1 hour	-
1000°C ~ 1300°C	10°C/min	30 min	-
1300°C ~ 1530°C	7°C/min	30 min	-
1530°C	-	2 hours	120 min
1,530°C ~ *R.T.	7°C/min	4 hours	-

*R.T. : Room Temperature

- Use the sintering schedule according to each system.
- Do not give any pressure when eliminating after sintering.
- Qualities can be improved if heat treatment is conducted at up to 1000°C prior to build-up after sintering has been completed.
- Use the build-up powder for zirconia use as recommended by the manufacturer.

E. Cautions

- (1) The user must understand the product thoroughly to prevent from any possible problems or disputes and the distributor must give an introduction regarding the procedure.
- (2) Keep an eye on the product for it is small-sized.
- (3) Avoid any contact from children and seniors.
- (4) Do not place it where it may become infected.
- (5) Do not give any pressure or impact.
- (6) It should be used only by dentists or specialists.
- (7) Recommended to follow the terms of the firing schedule table suggested.

F. Storage

- (1) Avoid high temperature and humidity and direct rays of the sun.
- (2) Keep it at normal temperature in a dry place and avoid any loss.

G. Transportation, Maintenance after use

- (1) Transport with cushioning material to buffer the impact.
- (2) Do not drop or throw.
- (3) The product is disposable and must not be reused.

[Pictograms]

	Manufacturer		Keep dry
	Do not reuse		Caution
	Date of Manufacture		Use by
	Batch Code		Catalogue Number
	Consult Instructions for Use		Non Sterile
	Do not use If package is damaged		Keep away from sunlight

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