

EVEREST Multilayer Universal Speed (UVS) combines our latest multilayer Zirconia technology to create a strong but transparent material with no compromise to shade stability in both speed and normal sintering. Suitable for single crowns and full arches, no variation in translucency or strength.

Translucency	Strength
48%	900MPa
47%	···· 1000MPa
46%	···· 1100MPa
45%	1150MPa

Shades A1, A2, A3, A3.5 B1, B2, B3 C1, C2, C3 D2, D3 BL1, BL3 Block SizeØ 95, 9812mm20mm14mm22mm16mm25mm18mm30mm

Comparison

Anterior

Crown



4 Hour

Why EVEREST Multilayer UVS?

- Fast sintering Zirconia
- 30 minute sintering cycle (based on single crown)
- Same results in both Normal and Speed sintering
- Same results in wide MAX temperature range (1500°C-1560°C)
- · Flexural strength, translucency, and shade stability
- For all indications (single crown to full arch, anterior and posterior



MULTILAYER UVS



EVEREST W ML UNIVERSAL SPEED Demension : Ø98

FDA ISO 13485:2016 9001:2015

Full Arch

U&C International Co., Ltd. 106-11, Beomjigi-ro, Danwon-gu Ansan-si, Gyeonggi-do, Korea www.uncint.com marketing@uncint.net

Indication

30 Minute

Single crown to full arch

(Anterior / Posterior)

INC

Posterior

Crown

One Day Zirconia

-EVECEST Volcan Speed

Paired with EVEREST Multilayer UVS, sinter a crown in only 30 minutes and full arch in only 4 hours.

- 30 minute speed sintering mode
- Normal sintering- 2 layer tray capacity (90 Ø)
- Speed sintering- Maximum 5ea crowns
- Shade stability and consistency
- High purity and heat durability





<u>everest</u>

Zirconia Pre-Polishing Felt

Paired with EVEREST Multilayer UVS and EVEREST Volcan Speed, finish restorations within one visit.

- Designed for polishing Zirconia before sintering
- Pre-polishing can save up to 30 minutes
- Takes only 30 seconds